




For Immediate Release: Nitinol Wire Reinforced Tubing

October 10, 2010

Contact: Ronald Lilly
Forefront Medical Technology
Farmington, CT 06034-1292
860.255.7610
rlilly@fmt.com.sg



Forefront Medical Technology has developed a new product range of specialty nitinol wire reinforced flexible tubing for high performance catheters and other unique delivery systems. Tubing constructions are available with medical grades of silicone, PVC or DEHP free PVC. Standard construction ranges are 5 mm to 20 mm outside diameters and up to 75 cm lengths. Custom constructions are available for unique application requirements.

They are a fully integrated contract manufacturer with global capabilities to meet customers' exacting manufacturing requirements. Forefront's expertise includes product design/development, rapid SLS prototyping, in-house tool making, injection molding, extrusion, class 10,000 clean room assembly, testing, packaging, sterilization and global logistics. In addition, they are ISO 9001:2008, ISO 13485:2003 and FDA CFR 21 Part 820 registered.

Forefront Medical's technical solutions reduce customers' time to market, ensure confidentiality of new designs and protection of intellectual property. For a confidential technical evaluation, please call 1 (860) 255-7610 (North America) / +65 6349 3888 (Asia) or on the web at <http://forefront.sg>.

#