



For Immediate Release: Specialty Silicone Tubing for Drug Delivery Pumps

February 18, 2010

Contact: Ronald Lilly
Forefront Medical Technology
Farmington, CT 06034-1292
860.255.7610
rlilly@fmt.com.sg



Forefront Medical Technology has developed a new specialty silicone tubing assembly for high performance drug delivery systems that mechanically locks the silicone pump tubing onto the connectors.

This innovative design eliminates the need for adhesives and the additional costs/concerns associated with adhesive bonding this critical assembly. Custom constructions are available for unique application requirements.

Forefront Medical is a fully integrated contract manufacturer with global capabilities to meet customers' exacting manufacturing requirements. Forefront's expertise focuses on injection molding, extrusion and blow molding, along with assembly and packaging. Complimenting their expertise are product design/development, rapid SLS prototyping, in-house tool making, class 100k molding and extrusion clean rooms, class 10k clean room assembly, product testing, sterilization and supply chain management. In addition, they are ISO 9001:2008, ISO 13485:2003 and FDA CFR 21 Part 820 registered.

Forefront Medical's technical solutions reduce customers' time to market, ensure confidentiality of new designs and protection of intellectual property. For a confidential technical evaluation, please call 1 (860) 255-7610 (North America) / +65 6349 3888 (Asia) or on the web at <http://forefront.sg>.

#

Xiamen, China

Singapore

Farmington, CT USA